STERILE A
Sterilized Using Aseptic Processing Techniques
Sterilizzato con Tecnica Asettica
Sterilisé à l'aide de techniques de traitement aseptique
Durch aseptische Herstellungsverfahren sterilisiert
Esterilizado con técnicas de procesamiento aséptico
Sterilized using aseptic processing methods
Стерилизовано с помощью асептической обработки
Esterilizado Utilizando Técnicas de Procesamiento Aséptico

Do not use if package is damaged
Non utilizzare se l'imballaggio è danneggiato
Ne pas utiliser si l'emballage est endommagé
Nicht verwenden, falls die Packung beschädigt ist
Nad utilizar se embalagem estiver danificada
Nie używać, jeśli opakowanie jest uszkodzone
Не использовать продукт с поврежденной упаковкой
No utilizar si el embalaje está dañado

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EN HYPURE™ OIL HEAVY

INTENDED USE

Hypure™ Oil is intended to be used as an overlay for culture of embryos, oocytes and sperm in ART and micromanipulation procedures.

COMPOSITION:

Pharmaceutical grade quality paraffin oil.

GENERAL INFORMATION:

Hypure™ Oil is pre-washed, ready to use mineral oil for use in IVF, ICSI and related artificial reproductive techniques.

The oil allows small volumes of medium to be used as it prevents the evaporation of media and maintain stable osmolality and pH in the culture medium.

INSTRUCTIONS FOR USE:

CAUTION: The detailed usage of Hypure™ Oil is determined based on each laboratory's procedure and protocols that are optimized for individual medical purposes.

Preparation :

- For culture, pre-equilibrate Hypure™ Oil as necessary, and the culture medium at 37°C in a CO₂ incubator for several hours (or overnight),
- For ICSI, warm Hypure™ Oil and HEPES buffered medium at 37°C.

Instructions for use :

1. Put culture medium on a dish or place micro drop of culture medium on a dish.
2. Gently pipette equilibrated Hypure™ Oil over the medium until they are completely covered.

NOTE : Washing of oil for embryo culture is not required.

QUALITY CONTROL SPECIFICATIONS:

The following tests were performed for each lot of this product:

- Sterility: SAL 10⁻³
- Endotoxin: < 0.1 EU/ml
- Mouse embryo assay: ≥ 80% blastocysts after 120 hours in culture
- Density: 0.86-0.88 g/ml
- Viscosity: 95-140 cP at 30°C
- Peroxide Value (POV): <0.1mEq/kg

Kitazato Corporation
100-10 Yanagishima, Fuji, Shizuoka 416-0932 Japan
TEL: +81-545-65-7122 FAX: +81-545-65-7128

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Copyright © KITAZATO CORPORATION.
All Rights Reserved.