
Vitrification Cryoprotocol for Oocyte and Embryo

Vitrification media

Code:VT601, VT602

Symbols on the Label

Ver.3

Sterilized using aseptic processing techniques	
STERILE A	Do not reuse
REF	Catalog number
LOT	Batch code
CE 2797	Use by
Manufacturer logo	Storage temperature
Do not re-sterilize	CE mark
Do not use if package is damaged	Manufacturer
Consult instructions for use	Do not re-sterilize
Authorized representative in the European Community	Do not use if package is damaged

EN: Intended use

This product is to be used for vitrification of oocytes (MI) and embryos.

Vitrification Media

No. 1 (Color code: white) Basic Solution (BS): 1x1.5 ml vial only for oocyte (MI) vitrification

No. 1 (Color code: blue) Equilibration Solution (ES): 1x1.5 ml vial only

No. 2 (Color code: green) Vitrification Solution (VS): 2x1.5 ml vials

Caution: Before use, check the specifications of the container and labeling (number marked on top of the cap, cap color, vial labeling color, solution name, and volume). If you notice anything unusual regarding the items mentioned above, do not use the product and please contact the distributor

Recommended equipment

Cryotop: 1 Cryotop stores up to 4 oocytes (MI) or 4 Embryos as a recommendation.
ReproPlate: with 6 wells
Cooling Rack

Caution: This product is sterilized. Please handle the product in a sterile field with a clean bench.

Instructions for use (IFU)**Preparation**

Fill 90% of the Cooling Rack with fresh liquid nitrogen.

Use 100 µl of sterilized saline with the thickness of zona pellucida and record it.

Bring BS, ES and VS to room temperature (23-27°C as recommended).

Use certified pipette as a handling tool, with a suitable internal diameter for oocytes or embryos.

The recommended internal diameters are as follows: 120 µm for oocytes (MI), 120 µm for pronuclear stage embryos, 150 to 180 µm for cleavage stage embryos, and 180 to 250 µm for blastocyst stage embryos.

Equilibration

Equilibration procedures for Oocyte (MI) and Embryo are different.

Equilibration of oocyte (MI)

Drop 20 µl of BS into first well and each 300 µl of VS into second and third well on the Repro Plate with pipette.

Step 0: Transfer the Embryo to the TOP of culture dish to the BOTTOM of BS.

Step 1: Immediately add ES 20 µl gently on the TOP of BS of the first well and leave it for 3 minutes.

Add another ES 20 µl gently on the TOP of the first well and leave it for 3 minutes.

Add another ES 240 µl gently on the TOP of the first well and leave it for 9 minutes.

Equilibration of embryos

Drop each 300 µl of ES into first well, VS into second and third well on the Repro Plate with pipette.

Step 0: Transfer the Embryo to the TOP of ES from ES of the culture dish. It will spontaneously begin to shrink and then gradually return to its original size by absorbing the ES solution (within 15 minutes).

Vitrification

Caution: The following steps from 1 to 9 should be completed between 60 and 90 seconds:

Aspirate the oocyte (MI)/embryo from ES with the tip of a pipette.

Transfer the oocyte (MI)/embryo to the TOP of center of the VS of second well.

Aspirate the oocyte (MI)/embryo with a pipette and blow it out. Repeat this process 3 times, changing the position in the VS of second well.

Transfer the oocyte (MI)/embryo to the VS of third well.

Change the position of the oocyte (MI)/embryo to the VS of third well with a pipette.

Place the oocyte (MI)/embryo by the black line on the Cryotop.

Make a planar drop.

Make sure the oocyte (MI)/embryo is on the Cryotop with a minimal volume of the VS and well (less than 0.1 µl) under a microscope.

Plug the Cryotop immediately in liquid nitrogen.

Store the Cryotop in a can and store it in a storage tank.

Quality control specification

The following tests were performed for each of this product:

Sterility by Sterility Test (EP)

Endotoxin by LL method.

Mouse Embryo Assay (One Cell)

pH (EP)

Osmolality (EP)

Storage instructions and stability

Store the vials at 2 to 8 °C.

This product is stable until the expiry date labeled on the vial.

Composition

-HEPES with Basic Culture Medium
-Ethylene Glycol
-Dimethylsulfoxide
-Trehalose
-Hydroxypropyl cellulose
-Gentamicin

Contraindications

Do not re-sterilize.

Do not use. Re-use may cause contamination.

Do not use solution that shows cloudiness or becomes discolored.

Do not use the product if you notice anything unusual regarding the solution on the label (number, color, name, volume).

Morphologically abnormal oocytes, embryos, or significantly poor grade oocytes/embryos are unsuitable for cryopreservation.

Observe all federal, state and local environmental regulations when discarding the product.

In case of use, dispose the product appropriately in a prescribed manner.

You shall be responsible for any problems caused by non-conformity to the present IFU.

Vitrification media contain an antibiotic gentamicin sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

NOTE: The long-term safety of cryopreservation technique and maximum storage in liquid nitrogen has not been established and unknown.

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