

Product Update

COMMUNICATION

Dear Distributors,

We are pleased to inform you that we are improving both the *Echo Marking* and the *Cork* of our Single Lumen OPU Needle range.

Echo Marking

It will be now smoother than the previous spiral carved type, and it will reach the tip of the needle, increasing the length of the current echogenic marking and extending the needle's overall echogenic surface. The new marking is identical to the one we implemented to the Double Lumen range.



Cork

The new *Cork* design allows the syringe to be directly connected to it, without the need of a cannula. This would make the procedure more hygienic and less time consuming. This improvement was applied into the Japanese market 3 years ago; not a single issue has been reported.



The improvement of the *Echo Marking* will be launched first, and then the new *Cork* design will follow. Our customer service department or your local distributor will inform you when your needle orders begin to include the new *Echo Marking*.

In the meantime, we remain at your service at orders@dibimed.com