

Hypure™ Oil: Heavy and Light



The best protection with the lowest endotoxicity.

Kitazato Hypure™ Oil is a high quality paraffin oil. It is **pre-washed with ultra-pure water** to reduce toxicity and remove endotoxins, therefore improving your clinical outcomes.

Our paraffin oils are indicated as **the best protective barrier to prevent the quick evaporation of the medium** and to keep specimens safe from external contaminants during in-vitro culture.

| Product | Code | Volume |
|-------------------|------------|--------|
| Hypure™ Oil Heavy | MOHVD-100P | 100 |
| | MOHVD-250P | 250 |
| | MOHVD-500P | 500 |
| Hypure™ Oil Light | MOLHD-100 | 100 |
| | MOLHD-100P | 100 |
| | MOLHD-250 | 250 |
| | MOLHD-500 | 500 |
| | MOLHD-500P | 500 |

KEY FEATURES

- Pharmaceutical grade quality **paraffin oil**.
- Allows small volumes of medium to be used as it prevents the evaporation of media.
- Maintains **stable osmolality** and pH in culture medium.
- Available in glass or PET bottles.
- **Hypure™ Oil Heavy** is pre-washed and offers extensive protection and insulation.
- **Hypure™ Oil Light** is double-washed and allows easy pipetting.

Available with different density and viscosity.